

Product Name: Recombinant Human ABCA9 Protein
Catalog #: PEH6749



Summary

Name	ABCA9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IUA7
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CFWFQHGRAN HVVLENETDS DPTPNDCFEP VSPEFCGKEA IRIKNLKKEY AGKCERVEAL KGVVFDIYEG QITALLGHSG AGKTTLLNIL SGLSVPTSGS VTVYNHTLSR MADIENISKF TGFCPQSNVQ FGFLT VKENL RLFAKIKGIL PHEVEKEVQR VVQEELEMENI QDILAQNLSG GQNRKLTFGI AILGDPQVLL LDEPTAGLDP LSRHRIWNLL KEGKSDRVIL FSTQFIDEAD ILADRKVFIS NGKLKAGSS LFLKKKWGIG YHLSLHLNER CDPESITSLV KQHISDAKLT AQSEEKLVYI LPLERTNKF ELYRDLDRCS NQGIEDYGVS ITTLNEVFLK LEGKSTIDES DIGIWGQLQT DGAKDIGSLV ELEQVLSSFH

image



Background

Alternative Names	ATP-binding cassette sub-family A member 9
--------------------------	--

Product Name: Recombinant Human ABCA9 Protein
Catalog #: PEH6749



Note

For Research Use Only , Not for Diagnostic Use