


Product Name: Recombinant Human ADCY6 Protein
Catalog #: PEH6753



Summary

Name	ADCY6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43306
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKEDMMFHKI YIQKHDNVS I LFADIEGFTS LASQCTAQEL VMTLNELFAR FDKLAENHC LRIKILGDCY YCVSGLPEAR ADHAHCCVEM GVDMIEAISL VREVTGVNVN MRVGIHSGRV HCGVLGLRKW QFDVWSNDVT LANHMEAGGR AGRIHITRAT LQYLNGDYEY EPGRGGERNA YLKEQHIETF LILGASQKRK EEKAMLAKLQ RTRANSM EGL MPRWVPDRAF SRTKDSKA FR QMGIDDSSKD NRGTQDALNP

image



Background

Alternative Names	Adenylate cyclase type 6 (EC 4.6.1.1) (ATP pyrophosphate-lyase 6) (Adenylate cyclase type VI) (Adenylyl cyclase 6) (Ca(2+)-inhibitable adenylyl cyclase)
--------------------------	--

Product Name: Recombinant Human ADCY6 Protein
Catalog #: PEH6753



Note

For Research Use Only , Not for Diagnostic Use