


Product Name: Recombinant Human CD7 Protein
Catalog #: PEH6787



Summary

Name	CD7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09564
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSPHCTTVPV GASVNITCST SGGLRGIYLR QLGPQPQDII YYEDGVVPTT DRRFRGRIDF SGSQDNLTIT MHRQLSDTG TYTCQAITEV NVYGSGLTVL VTEEQSQGWH RCSDAPPRAS ALPAPPTGSA LPDPQTASAL PDPPAASALP

image



Background

Alternative Names	T-cell antigen CD7 (GP40) (T-cell leukemia antigen) (T-cell surface antigen Leu-9) (TP41) (CD antigen CD7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use