


Product Name: Recombinant Human SLC25A12 Protein
Catalog #: PEH6791



Summary

Name	CMC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75746
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NSNPKIVQLL AGVADQTKDG LISYQEFLAF ESVLCAPDSM FIVAFQLFDK SGNGEVTFFEN VKEIFGQTII HHHIPFNWDC EFIRLHFGHN RKKHLNYTEF TQFLQELQLE HARQAFALKD KSKSGMISGL DFSDIMVTIR SHMLTPFVEE NLVSAAGGSI SHQVSFSYFN AFNSLLNNME LVRKIYSTLA GTRKDVEVTK EEFAQSAIRY GQVTPLEIDI LYQLADLYNA SGRLTLADIE RIAPLAEGAL PYNLAELQRQ QSPGLGRPIW

image



Background

Alternative Names	Calcium-binding mitochondrial carrier protein Aralar1 (Mitochondrial aspartate glutamate carrier 1) (Solute carrier family 25 member 12)
--------------------------	--

Product Name: Recombinant Human SLC25A12 Protein
Catalog #: PEH6791



Note

For Research Use Only , Not for Diagnostic Use