

Product Name: Recombinant Human DNAJC4 Protein
Catalog #: PEH6806



Summary

Name	DNJC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NNZ3
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPSTYYELLG VHPGASTEEV KRAFFSKSKE LHPDRDPGPNP SLHSRFVELS EAYRVLREQ SRRSYDDQLR SGSPPKSPRT TVHDKSAHQ T

image



Background

Alternative Names DnaJ homolog subfamily C member 4 (DnaJ-like protein HSPF2) (Multiple endocrine neoplasia type 1 candidate protein number 18)

Note

For Research Use Only , Not for Diagnostic Use