


Product Name: Recombinant Human CA14 Protein
Catalog #: PEH6817



Summary

Name	CAH14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9ULX7
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HWTYEGPHGQ DHWPASYPEC GNNAQSPIDI QTDSVTDFPD LPALQPHGYD QPGTEPLDLH NNGHTVQLSL PSTLYLGGLP RKYVAAQLHL HWGQKGSPGG SEHQINSEAT FAELHIVHYD SDSYDSLSEA AERPQGLAVL GILIEVGETK NIAYEHILSH LHEVRHKDQK TSVPPFNLRE LLPKQLGQYF RYNGSLTTPP CYQSVLWTVF YRRSQISMEQ LEKLQGTLS TEEPSKLLV QNYRALQLN QRMVFASFIQ AGSSYTTGEM

image



Background

Alternative Names	Carbonic anhydrase 14 (EC 4.2.1.1) (Carbonate dehydratase XIV) (Carbonic anhydrase XIV) (CA-XIV)
--------------------------	--

Product Name: Recombinant Human CA14 Protein
Catalog #: PEH6817



Note

For Research Use Only , Not for Diagnostic Use