


Product Name: Recombinant Human ALDH3B1 Protein
Catalog #: PEH6831



Summary

Name	AL3B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P43353
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKHLTPVTLE LGGKNPCYVD DNCDPQTVAN RVAWFRYFNA GQTCVAPDYV LCSPQMQLERL LPALQSTITR FYGDDPQSSP NLGRIINQKQ FQRLRALLGC GRVAIGGQSD ESDRYIAPTV LVDVQEMEPV MQEEIFGPIL PIVNVQSLDE AIEFINRREK PLALYAFSNS SQVVKRVLQ TSSGGFCGND GFMHMTLASL PFGGVGASGM GRYHGKFSFD TFSHHRACLL RSPGMEKLNA LRYPPQSPRR

image



Background

Alternative Names	Aldehyde dehydrogenase family 3 member B1 (EC 1.2.1.28) (EC 1.2.1.5) (EC 1.2.1.7) (Aldehyde dehydrogenase 7)
--------------------------	--

Note

Product Name: Recombinant Human ALDH3B1 Protein
Catalog #: PEH6831



For Research Use Only , Not for Diagnostic Use