


Product Name: Recombinant Human ALDH3B1 Protein
Catalog #: PEH6832



Summary

Name	AL3B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P43353
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKHLTPVTLE LGGKNPCYVD DNCDPQTVAN RVAWFRYFNA GQTCVAPDYV LCSPERMQERL LPALQSTITR FYGDDPQSSP NLGRIINQKQ FQRLRALLGC GRVAIGGQSD ESDRYIAPT V LVDVQEMEPV MQEEIFGPIL PIVNVQSLDE AIEFINRREK PLALYAFSNS SQVVKRVL TQ TSSGGFCGND GFMHMTLASL PFGGVGASGM GRYHGKFSFD TFSHHRACLL RSPGMEKLNA LRYPPQSPRR

image



Background

Alternative Names Aldehyde dehydrogenase family 3 member B1 (EC 1.2.1.28) (EC 1.2.1.5) (EC 1.2.1.7) (Aldehyde dehydrogenase 7)

Note

Product Name: Recombinant Human ALDH3B1 Protein
Catalog #: PEH6832



For Research Use Only , Not for Diagnostic Use