


Product Name: Recombinant Human AKR7A2 Protein
Catalog #: PEH6836



Summary

Name	ARK72_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43488
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MLSAASRVVS RAAVHICALRS PPPEARALAM SRPPPPRVAS VLGTMEMGRR MDAPASAAAV RAFLERGHTE LDATFMYSDG QSETILGGLG LGLGGGDCRV KIATKANPWD GKSLKPDSVR SQLETSLKRL QCPQVDLFYL HAPDHGTPVE ETLHACQRLH QEGKFVELGL SNYASWEVAE ICTLCKSNGW ILPTVYQGMV NATTRQVETE LFPCLRHFGL RFYAYNPLAG GLLTGKYKYE DKDGKQPVGR FFGNSWAETY RNRFWKEHHF EAIALVEKAL QAAYGASAPS VTSAAALRWMY HHSQLQGAHG DAVILGMSSL EQLEQNLAAT EEGPLEPAVV DAFNQAWHLV AHECPNYFR

image



Background

Alternative Names	Aflatoxin B1 aldehyde reductase member 2 (EC 1.1.1.n11) (AFB1 aldehyde
--------------------------	--

Product Name: Recombinant Human AKR7A2 Protein
Catalog #: PEH6836



reductase 1) (AFB1-AR 1) (Aldoketoreductase 7) (Succinic semialdehyde reductase) (SSA reductase)

Note

For Research Use Only , Not for Diagnostic Use