


Summary

Name	MDHC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40925
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSEPIRVLVT GAAGQIAYSL LYSIGNSVF GKDQPIILVL LDITPMMGVL DGVLMELQDC ALPLLKDZIA TDKEDVAFKD LDVAILVGSM PRREGMERKD LLKANVKIFK SQGAALDKYA KKSVMVIVVG NPANTNCLTA SKSAPSIPKE NFSCLTRLDH NRAKAQIALK LGVTANDVKN VIIWGNHSST QYPDVNHAKV KLOGKEVGVY EALKDDSWLK GEFVTTVQQR GAAVIKARKL SSAMSAAKAI CDHVRDIWFG TPEGEFVSMG VISDGNSYGV PDDLVSFPV VIKNKTWKVF EGLPINDFSR EKMDLTAKEL TEEKESAFEF LSSA

image



Background

Alternative Names	Malate dehydrogenase, cytoplasmic (EC 1.1.1.37) (Cytosolic malate dehydrogenase) (Diodophenylpyruvate reductase) (EC 1.1.1.96)
--------------------------	--

Product Name: Recombinant Human MDH1 Protein
Catalog #: PEH6837



Note

For Research Use Only , Not for Diagnostic Use