


Product Name: Recombinant Human RPS7 Protein
Catalog #: PEH6838



Summary

Name	RS7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62081
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MFSSSAKIVK PNGEKPDEFE SGISQALLEL EMNSDLKAQL RELNITAAKE IEVGGGRKAI IIFVPVQLK SFQKIQVRLV RELEKKFSGK HVVFIAQRRI LPKPTRKSRT KNKQKRPRSR TLTAVHDAIL EDLVFPSEIV GKRIKRVKLDG SRLIKVHLDK AQQNNVEHKV ETFSGVYKKL TGKDVNFEFP EFQL

image



Background

Alternative Names 40S ribosomal protein S7 (Small ribosomal subunit protein eS7)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RPS7 Protein
Catalog #: PEH6838

