


Product Name: Recombinant Human VTA1 Protein
Catalog #: PEH6841



Summary

Name	VTA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NP79
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAALAPLPL PAQFKSIQHH LRTAQEHDKR DPVVAYYCRL YAMQTGMKID SKTPECRKFL SKLMDQLEAL KKQLGDNEAI TQEIVGCAHL ENYALKMFLY ADNEDRAGRF HKNMIKSFYT ASLLIDVITV FGELTDENVK HRKYARWKAT YIHNCLKNGE TPQAGPVGIE EDNDIEENED AGAASLPTQP TQPSSSSTYD PSNMPSGNYT GIQIPGAHA PANTPAEVPH STGVASNTIQ PTPQTIPAI PALFNTISQG DVRLTPEDFA RAQKYCKYAG SALQYEDVST AVQNLQKALK LLTTGRE

image



Background

Alternative Names Vacuolar protein sorting-associated protein VTA1 homolog (Dopamine-responsive gene 1 protein) (DRG-1) (LYST-interacting protein 5) (LIP5) (SKD1-

Product Name: Recombinant Human VTA1 Protein
Catalog #: PEH6841



binding protein 1) (SBP1)

Note

For Research Use Only , Not for Diagnostic Use