


Product Name: Recombinant Human AS3MT Protein
Catalog #: PEH6844



Summary

Name	AS3MT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HBK9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CQVIYNGGIT GHEKELMFDA NTFKEGEIV EVDEETAAIL KNSRFAQDFL IRPIGEKLPT SGGCSALELK DIITDPFKLA

image



Background

Alternative Names Arsenite methyltransferase (EC 2.1.1.137) (Methylarsonite methyltransferase) (S-adenosyl-L-methionine:arsenic(III) methyltransferase)

Note

For Research Use Only , Not for Diagnostic Use