


Product Name: Recombinant Human CYP11B1 Protein
Catalog #: PEH6845



Summary

Name	C11B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	P15538
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTRAARVPRT VLPFEAMPRR PGNRWLRLQ IWREQGYEDL HLEVHQTFQE LGPIFRYDLG GAGMVCVMLP EDVEKLQQVD SLPHRMSLE PWVAYRQHRG HKCGVFLNG PEWRFNRLRL NPEVLSPNAV QRFLPMVDAV ARDFSQALKK KVLQNARGSL TLDVQPSIFH YTIEASNLAL FGERLGLVGH SPSSASLNFL HALEVMFKST VQLMFMPRSL SRWTSPKVWK EHFEAWDCIF QYGDNCIQKI YQELAFSRPQ QYTSIVAELL LNAELSPDAI KANSMELTAG SVDTTVFPLL MTLFELARNP NVQQALRQES LAAAASISEH PQKATTELPL LRAALKETLR LYPVGLFLER VASSDLVLQN YHIPAGTLVR VFLYSLGRNP ALFPRPERYN PQRWLDIRGS GRNFYHVPFG FGMRQCLGRR LAEAEMLLLL HHVLKHLQVE TLTQEDIKMV YSFILRPSMF PLLTFRAIN

image



Background

Product Name: Recombinant Human CYP11B1 Protein
Catalog #: PEH6845



Alternative Names

Cytochrome P450 11B1, mitochondrial (CYPXIB1) (Cytochrome P-450c11) (Cytochrome P450C11) (Steroid 11-beta-hydroxylase, CYP11B1) (EC 1.14.15.4)

Note

For Research Use Only , Not for Diagnostic Use