


Product Name: Recombinant Human FAM166B Protein
Catalog #: PEH6846



Summary

Name	F166B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	A8MTA8
Host	E.coli
Species	Human
Molecular Mass	70 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GYTGHCP LLR FSVGQTYGQV TGQLLRGPPG LAWPPVHRTL LPPIRPPRSP EVPRESLPVR RGQERLSSSM IPGYTGFVPR AQFIFAKNCS QVWAEALSDF THLHEKQGSE ELPKEAKGRK DTEKDQVPEP EGQLEEPTLE VVEQASPYSM DDRDP RKFFM SGFTGYVPCA RFLFGSSFPV LTNQALQEFQ QKHSPGSAQD PKHLPPLPRT YPQNLGLLPN YGGYVPGYKF QFGHTFGHLT

image



Background

Alternative Names Protein FAM166B

Note

Product Name: Recombinant Human FAM166B Protein
Catalog #: PEH6846



For Research Use Only , Not for Diagnostic Use