


Product Name: Recombinant Human ID4 Protein
Catalog #: PEH6847



Summary

Name	ID4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P47928
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKAVSPVRPS GRKAPSGCGG GELALRCLAE HGHSLGGSAA AAAAAAARC KAAEAAADEP ALCLQCDMND CYSRLRRLVP TIPPNNKVKSK VEILQHVIDY ILDQLALET HPALLRQPPP PAPPHPAGT CPAAPPRTPL TALNTDPAGA V NKQGDSILC R

image



Background

Alternative Names DNA-binding protein inhibitor ID-4 (Class B basic helix-loop-helix protein 27) (bHLHb27) (Inhibitor of DNA binding 4) (Inhibitor of differentiation 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ID4 Protein
Catalog #: PEH6847

