


Product Name: Recombinant Human COPE Protein
Catalog #: PEH6848



Summary

Name	COPE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14579
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAPPAPGPAS GGSGEVDRLF DVKNFYIGS YQQCINEAQR VKLSSPERDV ERDVFLYRAY LAQRKFGVVL DEIKPSSAPE LQAVRMFADY LAHESRRDSI VAELDREMSR SVDVTNTTFL LMAASIYLDH QNPDAALRAL HQGDSLECTA MTVQILLKLD RLDLARKELK RMQDLDEDAT LTQLATAWVS LATGGEKLQD AYYIFQEMAD KCSPTLLLLN GQAACHMAQG RWEAAEGLLQ EALDKDSGYP ETLVNLIVLS QHLGKPPEVT NRYLSQLKDA HRSHPIKEY QAKENDFDRL VLQYAPSA

image



Background

Alternative Names	Coatmer subunit epsilon (Epsilon-coat protein) (Epsilon-COP)
--------------------------	--

Product Name: Recombinant Human COPE Protein
Catalog #: PEH6848



Note

For Research Use Only , Not for Diagnostic Use