


Product Name: Recombinant Human YOD1 Protein
Catalog #: PEH6851



Summary

Name	OTU1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VVQ6
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MFGPAKGRHF GVHPAPGFPG GVSQQAAGTK AGPAGAWPVG SRTDTMWRLR CKAKDGTHTL QGLSSRTRVR ELQQIAAIT GIAPGGQRIL VGYPPECLDL SNGDTILEDL PIQSGDMLII EEDQTRPRSS PAFTKRGASS YVRETLPLVT RTVVPADNSC LFTSVYYVVE GGVLNPACAP EMRRLIAQIV ASDPDFYSEA ILGKTNQEYC DWIKRDDTWG GAIEISILSK FYQCEICVVD TQTVRIDRFG EDAGYTKRVL LIYDGIHYDP LQRNFPDPDT PPLTIFSSND DIVLVQALEL ADEARRRRQF TDVNRFTLRC MVCQKGLTGQ AEAREHAKET GHTNFGEV

image



Background

Alternative Names Ubiquitin thioesterase OTU1 (EC 3.4.19.12) (DUBA-8) (HIV-1-induced protease 7) (HIN-7) (HsHIN7) (OTU domain-containing protein 2)

Product Name: Recombinant Human YOD1 Protein
Catalog #: PEH6851



Note

For Research Use Only , Not for Diagnostic Use