


**Product Name: Recombinant Human DCPS Protein**  
**Catalog #: PEH6852**



## Summary

<b>Name</b>	DCPS_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96C86
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MADAAPQLGK RKRELDVEEA HAASTEEKEA GVGNGTCAPV RLPFSGFRLQ KVLRESARDK IIFLHGKVNE ASGDGDGEDA VVILEKTPFQ VEQVAQLLTG SPELQLQFSN DIYSTYHLFP PRQLNDVKTT VVYPATEKHL QKYLRQDLRL IRETGDDYRN ITLPHLESQS LSIQWVYNIL DKKAEADRIV FENPDPSDGF VLIPDLKWNQ QQLDDLYLIA ICHRRGIRSL RDLTPEHLPL LRNILHQQQE AILQRYRMKG DHLRVYLHYL PSYYHLHVHF TALGFEAPGS GVERAHLAE VIENLECDPR HYQQRILTFA LRADDPLLKL LQEAQQS

## image



## Background

**Alternative Names** m7GpppX diphosphatase (EC 3.6.1.59) (DCS-1) (Decapping scavenger enzyme) (Hint-related 7meGMP-directed hydrolase) (Histidine triad

**Product Name: Recombinant Human DCPS Protein**  
**Catalog #: PEH6852**



---

nucleotide-binding protein 5) (Histidine triad protein member 5) (HINT-5)  
(Scavenger mRNA-decapping enzyme Dcp5)

**Note**

For Research Use Only , Not for Diagnostic Use