

**Product Name: Recombinant Human PDLIM1 Protein**  
**Catalog #: PEH6855**



## Summary

<b>Name</b>	PDLI1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O00151
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MTTQQIDLQG PGPWGFRLVG GKDFEQPLAI SRVTPGSKAA LANLCIGDVI TAIDGENTSN MTHLEAQNRI KGCTDNLTLT VARSEHKVWS PLVTEEGKRH PYKMNLASEP QEVLHIGSAH NRSAMPFTAS PASSTTARVI TNQYNNPAGL YSSNISNFN NALESKTAAS GVEANSRPLD HAQPPSSLVI DKESEVYKML QEKQELNEPP KQSTSFLVLQ EILESEKGD PNKPSGFRSV KAPVTKVAAS IGNAQKLPMC DKCGTGIVGV FVKLDRHRH PECYVCTDCG TNLKQKGHFF VEDQIYCEKH ARERVTPPEG YEVVTVFPK

## image



## Background

**Alternative Names** PDZ and LIM domain protein 1 (C-terminal LIM domain protein 1) (Elfin) (LIM domain protein CLP-36)

**Product Name: Recombinant Human PDLIM1 Protein**  
**Catalog #: PEH6855**

---



**Note**

For Research Use Only , Not for Diagnostic Use