


Product Name: Recombinant Human TPM3 Protein
Catalog #: PEH6861



Summary

Name	TPM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P06753
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MMEAIIKKMQ MLKLDKENAL DRAEQAEAEQ KQAEERSKQL EDELAAMQKK LKGTEDELDK YSEALKDAQE KLELAEKKAA DAEAEVASLN RRIQLVEEEL DRAQERLATA LQKLEEAeka ADESERGMKV IENRALKDEE KMELQEIQLK EAKHIAEEAD RKYEEVARKL VIIEGDLERT EERAELAESK CSELEEEELKN VTNNLKSLEA QAEKYSQKED KYEEEIKILT DKLKEAETRA EFAERSVAKL EKTIDDELE LYAQLKLYKA ISEELDHAN DMTSI

image



Background

Alternative Names Tropomyosin alpha-3 chain (Gamma-tropomyosin) (Tropomyosin-3) (Tropomyosin-5) (hTM5)

Product Name: Recombinant Human TPM3 Protein
Catalog #: PEH6861



Note

For Research Use Only , Not for Diagnostic Use