


Product Name: Recombinant Human PCMT1 Protein
Catalog #: PEH6867



Summary

Name	PIMT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P22061
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAWKSGGASH SELIHNLRKN GIIKTDKVFE VMLATDRSHY AKCNPYMDSP QSIGFQATIS APHMHAYALE LLFDQLHEGA KALDVGSGSG ILTACFARMV GCTGKIVIGID HIKELVDDSV NNVRKDDPTL LSSGRVQLVV GDGRMGYAE APYDAIHVGA AAPVVPQALI DQLKPGGRLI LPVGPAGGNQ MLEQYDKLQD GSIKMKPLMG VIYVPLTDKE KQWSRWK

image



Background

Alternative Names	Protein-L-isoaspartate(D-aspartate) O-methyltransferase (PIMT) (EC 2.1.1.77) (L-isoaspartyl protein carboxyl methyltransferase) (Protein L-isoaspartyl/D-aspartyl methyltransferase) (Protein-beta-aspartate methyltransferase)
--------------------------	---

Product Name: Recombinant Human PCMT1 Protein
Catalog #: PEH6867



Note

For Research Use Only , Not for Diagnostic Use