


Product Name: Recombinant Human MTHFS Protein
Catalog #: PEH6873



Summary

Name	MTHFS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49914
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAAVSSAK RSLRGELKQR LRAMSAEERL RQSRVLSQKV IAHSEYQSKS RISIFLSMQD EIETEEIHKD IFQRGKICFI PRYRFQSNHM DMVRIESPEE ISLLPKTSWN IPQPGEGDVR EEALSTGGLD LIFMPGLGFD KHGNRLGRGK GYDAYLKRC LQHQEVKPYT LALAFKEQIC LQVPVNENDM KVDEVLYEDS STA

image



Background

Alternative Names	5-formyltetrahydrofolate cyclo-ligase (EC 6.3.3.2) (5,10-methenyl-tetrahydrofolate synthetase) (MTHFS) (Methenyl-THF synthetase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MTHFS Protein
Catalog #: PEH6873

