


Product Name: Recombinant Human GATM Protein
Catalog #: PEH6875



Summary

Name	GATM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P50440
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PKDCPVSSYN EWDPLEEVIV GRAENACVPP FTIEVKANTY EKYWPFYQKQ GGHYFPKDHL KKAVAEIEEM CNILKTEGVT VRRPDPIDWS LKYKTPDFES TGLYSAMPRD ILIVVGNEII EAPMAWRSRF FEYRAYRSII KDYFHRGAKW TTAPKPTMAD ELYNQDYPIH SVEDRHKLAA

image



Background

Alternative Names Glycine amidinotransferase, mitochondrial (EC 2.1.4.1) (L-arginine:glycine amidinotransferase) (Transamidinase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GATM Protein
Catalog #: PEH6875

