


Product Name: Recombinant Human BCAT1 Protein
Catalog #: PEH6876



Summary

Name	BCAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P54687
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKDCSNGCSA ECTGEGGSKE VVGTFKAKDL IVTPATILKE KPDPNNLVFG TVFTDHMLTV EWSSEFGWEK PHIKPLQNL LHPGSSALHY AVELFEGLKA FRGVDNKIRL FQPNLNMDRM YRSAVRATLP VFDKEELLEC IQQLVKLDQE WVPYSTSASL YIRPTFIGTE PSLGVKKPTK ALLFVLLSPV GPYFSSGTFN PVSLWANPKY VRAWKGGTGD CKMGGNYGSS LFAQCEAVDN GCQQVLWLYG EDHQITEVGT MNLFYWINE DGEELATPP LDGIILPGVT RRCILDLAHQ WGEFKVSERY LTMDDLTTAL EGNRVREMFG SGTACVVCV SDILYKGETI HIPTMENGPK LASRILSKLT DIQYGREESD WTIVLS

image



Background

Alternative Names Branched-chain-amino-acid aminotransferase, cytosolic (BCAT(c)) (EC

Product Name: Recombinant Human BCAT1 Protein
Catalog #: PEH6876



2.6.1.42) (Protein ECA39)

Note

For Research Use Only , Not for Diagnostic Use