


Product Name: Recombinant Human SCRN1 Protein
Catalog #: PEH6879



Summary

Name	SCRN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12765
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAAPPSYCF VAFPPRAKDG LVVFGKNSAR PRDEVQEVVY FSAADHEPES KVECTYISID QVPRTYAIMI SRPAWLWGAE MGANEHGVC I ANEAINTREP AAEIEALLGM DLVRLGLERG ETAKEALDVI VSLLEEHGQG GNYFEDANSC HSFQSAYLIV DRDEAWVLET IGKYWAAEKV TEGVRCICSQ LSLTTKMDAE HPELRSYAQS QGWWTGEGEF NFSEVFSPVE DHLDCGAGKD SLEKQEE SIT VQTM MNTLRD KASGVCIDSE FFLTASGVS VLPQNRSSPC IHYFTGTPDP SRSIFKPFIF VDDVKLVPKT QSPCFGDDDP AKKEPRFQEK PDRRHELYKA HEWARAIIES DQEQRKRLRS TMLELEKQGL EAMEEILTSS EPLDPAEVDG LFYDCVDTEI KFFK

image



Background

Product Name: Recombinant Human SCRN1 Protein
Catalog #: PEH6879



Alternative Names

Secernin-1

Note

For Research Use Only , Not for Diagnostic Use