


Product Name: Recombinant Human SEPTIN2 Protein
Catalog #: PEH6882



Summary

Name	SEPT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q15019
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERERLKKRIL DEIEHNIKI YHLPDAESDE DEDFKEQTRL LKASIPFSVV GSNQLIEAKG KKVGRGLYPW GVVEVENPEH NDFLKLRTML ITHMQDLQEV TQDLHYENFR SERLKRGGGRK VENEDMNKDQ ILLEKEAELR RMQEMIARMQ AQMQMQMGGG DGDGGALGHH

image



Background

Alternative Names	Septin-2 (Neural precursor cell expressed developmentally down-regulated protein 5) (NEDD-5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SEPTIN2 Protein
Catalog #: PEH6882

