

Product Name: Recombinant Human PTGR2 Protein
Catalog #: PEH6885




Summary

Name	PTGR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N8N7
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MIVQRVVLNS RPGKNGNPVA ENFRMEEVYL PDNINEGQVQ VRTLKLSVDP YMRCRMNEDT GTDYITPWQL SQVVDGGGIG IIESKHTNL TKGDFVTSFY WPWQTKVILD GNSLEKVPDQ LVDGHLSYFL GAIGMPGLTS LIGIQEKGHI TAGSNKTMVV SGAAGACGSV AGQIGHFLGC SRVVGICGTH EKCILLTSEL GFDAAINYKK DNVAEQLRES CPAGVDVYFD NVGGNISDTV ISQMNENSHI ILCGQISQYN KDVPYPPPLS PAIEAIQKER NITRERFLVL NYKDKFEPGI LQLSQWFKEG KLKIKETVIN GLENMGAAFQ SMMTGGNIGK QIVCISEEIS L

image



SDS



WB

Background

Alternative Names	Prostaglandin reductase 2 (PRG-2) (EC 1.3.1.48) (15-oxoprostaglandin 13-reductase) (Zinc-binding alcohol dehydrogenase domain-containing protein)
--------------------------	---

Product Name: Recombinant Human PTGR2 Protein
Catalog #: PEH6885



1)

Note

For Research Use Only , Not for Diagnostic Use