


Summary

Name	BHMT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q93088
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPPVGGKKAK KGILERLNAG EIVIGDGGFV FALEKRGYVK AGPWTPAAV EHPEAVRQLH REFLRAGSNV MQTFTFYASE DKLENRGNVY LEKISGQEVN EAACDIARQV ADEGDALVAG GVSQTPSYLS CKSETEVKKV FLQQLEVFMK KNVDFLIAEY FEHVEEAVWA VETLIASGKP VAATMCIGPE GDLHGVPPE CAVRLVKAGA SIIGVNCHFD PTISLKTVKL MKEGLEAARL KAHLMSQPLA YHTPDCNKQG FIDLPEFPFG LEPRVATRWD IQKYAREAYN LGVRYIGGCC GFEPYHIRAI AEELAPERGF LPPASEKHGS WGSGLDMHTK PWVRARARKE YWENLRIASG RPYNPSMSKP DGWGVTKGTA ELMQQKEATT EQQLKELFEK QKFKSQ

image



Background

Product Name: Recombinant Human BHMT Protein
Catalog #: PEH6886



Alternative Names

Betaine--homocysteine S-methyltransferase 1 (EC 2.1.1.5)

Note

For Research Use Only , Not for Diagnostic Use