



**Product Name: Recombinant Human PSMB7 Protein**  
**Catalog #: PEH6887**



## Summary

<b>Name</b>	PSB7_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q99436
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	26 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	TTIAGVVYKD GIVLGADTRA TEGMVVADKN CSKIHFISPN IYCCGAGTAA DTDMTTQLIS SNLELHSLST GRLPRVVTAN RMLKQMLFRY QGYIGAALVL GGVDVTGPHL YSIYPHGSTD KLPYVTMGSG SLAAMAVFED KFRPDMEEEE AKNLVSEAIA AGIFNDLGSG SNIDLCVISK NKLDFLRPYT VPNKKGTRLG RYRCEKGTTA VLTEKITPLE IEVLEETVQT MDTS

## image



## Background

**Alternative Names** Proteasome subunit beta type-7 (EC 3.4.25.1) (Macropain chain Z)  
(Multicatalytic endopeptidase complex chain Z) (Proteasome subunit Z)

## Note

**Product Name: Recombinant Human PSMB7 Protein**  
**Catalog #: PEH6887**



---

For Research Use Only , Not for Diagnostic Use