


Product Name: Recombinant Human COPS8 Protein
Catalog #: PEH6890



Summary

Name	CSN8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q99627
Host	E.coli
Species	Human
Molecular Mass	62 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FSFKLLDQC ENQELEAPGG IATPPVYGQL LALYLLHNDM NNARYLWKRI PPAIKSANSE LGGIWSVGQR IWQRDFPGIY TTINAHQWSE TVQPIMEALR DATRRRAFAL VSQAYTSIIA DDFAAFVGLP VEEAVKGILE QGWQADSTTR MVLPRKPVAG ALDVSNKFI

image



Background

Alternative Names COP9 signalosome complex subunit 8 (SGN8) (Signalosome subunit 8) (COP9 homolog) (hCOP9) (JAB1-containing signalosome subunit 8)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human COPS8 Protein
Catalog #: PEH6890

