


Product Name: Recombinant Human COPS8 Protein
Catalog #: PEH6891



Summary

Name	CSN8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q99627
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPVAVMAESA FSFKLLDQC ENQELEAPGG IATPPVYGQL LALYLLHNDM NNARYLWKRI PPAIKSANSE LGGIWSVGQR IWQRDFPGIY TTINAHQWSE TVQPIMEALR DATRRRAFAL VSQAYTSIIA DDFAAFVGLP VEEAVKGILE QGWQADSTTR MVLPRKPVAG ALDVSNKFI PLSEPAPVPP IPNEQQLARL TDYVAFLEN

image



Background

Alternative Names COP9 signalosome complex subunit 8 (SGN8) (Signalosome subunit 8) (COP9 homolog) (hCOP9) (JAB1-containing signalosome subunit 8)

Note

Product Name: Recombinant Human COPS8 Protein
Catalog #: PEH6891



For Research Use Only , Not for Diagnostic Use