


Product Name: Recombinant Human NPL Protein
Catalog #: PEH6893



Summary

Name	NPL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXD5
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAFPKKKLQG LVAATITPMT ENGEINFSVI GQYVDYLVKE QGVKNIFVNG TTGEGLSLSV SERRQVAEEW VTKGKDKLDQ VIIHVGALS L KESQELAQHA AEIGADGIAV IAPFFLKPWT KDILINFLKE VAAAAPALPF YYYHIPALTG VKIRAEELLD GILDKIPTFQ GLKFSDTDLL DFGQCVDQNR QQQFAFLFGV DEQLLSALVM GATGAVGSTY NYLGKKTNQM LEAFEQKDFS LALNYQFCIQ RFINFVVKLG FGVSQTKAIM TLVSGIPMGP PRLPLQKASR EFTDSAEAKL KSLDFLSFTD LKDG NLEAGS

image



Background

Alternative Names	N-acetylneuraminate lyase (NALase) (EC 4.1.3.3) (N-acetylneuraminate pyruvate-lyase) (N-acetylneuraminic acid aldolase) (Sialate lyase) (Sialate-
--------------------------	---

Product Name: Recombinant Human NPL Protein
Catalog #: PEH6893



pyruvate lyase) (Sialic acid aldolase) (Sialic acid lyase)

Note

For Research Use Only , Not for Diagnostic Use