


Product Name: Recombinant Human GLTP Protein
Catalog #: PEH6896



Summary

| | |
|------------------------------|--|
| Name | GLTP_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9NZD2 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 25 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | MALLAEHLLK PLPADKQIET GPFLEAVSHL PPFDCLGSP VFTPIKADIS GNITKIKAVY DTNPAKFRTL QNILEVEKEM YGAEWPKVGA TLALMWLKRK LRFIQVFLQS ICDGERDENH PNLIRVNATK AYEMALKKYH GWIVQKIFQA ALYAAPYKSD FLKALSKGQN VTEEECLEKI RLFLVNYTAT IDVIYEMYTQ MNAELNYKV |

image



Background

Alternative Names Glycolipid transfer protein (GLTP)

Note

Product Name: Recombinant Human GLTP Protein
Catalog #: PEH6896



For Research Use Only , Not for Diagnostic Use