


Product Name: Recombinant Human ARPC1B Protein
Catalog #: PEH6899



Summary

Name	ARC1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O15143
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQIAICPNNH EVHIYEKSGA KWTKVHELKE HNGQVTGIDW APESNRIVTC GTDRNAYVWT LKGRTWKPTL VILRINRAAR CVRWAPNENK FAVGSGSRVI SICYFEQEND WWVCKHIKKP IRSTVLSLDW HPNNVLLAAG SCDFKCRIFS AYIKEVEERP APTPWGSKMP FGELMFESS SCGWVHGVCF SASGSRVAWV SHDSTVCLAD ADKKMAVATL ASETLPLLAL TFITDNSLVA AGHDCFPVLF TYDAAAGMLS FGGRLDVPKQ SSQRGLTARE RFQNLDDKKAS SEGGAAGAG LDSLHKNSVS

image



Background

Alternative Names	Actin-related protein 2/3 complex subunit 1B (Arp2/3 complex 41 kDa subunit) (p41-ARC)
--------------------------	--

Product Name: Recombinant Human ARPC1B Protein
Catalog #: PEH6899



Note

For Research Use Only , Not for Diagnostic Use