


Product Name: Recombinant Human ARL6IP5 Protein
Catalog #: PEH6904



Summary

Name	PRAF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75915
Host	E.coli
Species	Human
Molecular Mass	7 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRNLKNKLEN KMEGIGLKRT PMGIVLDALE QQEEGINRLT

image



Background

Alternative Names	PRA1 family protein 3 (ADP-ribosylation factor-like protein 6-interacting protein 5) (ARL-6-interacting protein 5) (Aip-5) (Cytoskeleton-related vitamin A-responsive protein) (Dermal papilla-derived protein 11) (GTRAP3-18) (Glutamate transporter EAAC1-interacting protein) (JM5) (Prenylated Rab acceptor protein 2) (Protein JWa) (Putative MAPK-activating protein PM27)
--------------------------	--

Note

Product Name: Recombinant Human ARL6IP5 Protein
Catalog #: PEH6904



For Research Use Only , Not for Diagnostic Use