


Product Name: Recombinant Human HBD Protein
Catalog #: PEH6908



Summary

Name	HBD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P02042
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVNALWGKVN VDAVGGEALG RLLVVYPWTQ RFFESFGDLS SPDAVMGNPK VKAHGKKVLG AFSDGLAHLN NLKGTFSQLS ELHCDKLHVD PENFRLLGNV LVCVLARNFG KEFTPMQAA YQKVVAGVAN

image



Background

Alternative Names Hemoglobin subunit delta (Delta-globin) (Hemoglobin delta chain)

Note

For Research Use Only , Not for Diagnostic Use