


Product Name: Recombinant Human RPLP2 Protein
Catalog #: PEH6909



Summary

Name	RLA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05387
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MRYVASYLLA ALGGNSSPSA KDIKKILDSV GIEADDDRLN KVISLNGKN IEDVIAQGIG KLASVPAGGA VAVSAAPGSA APAAGSAPAA AEEKKDEKKE ESESDDDMG FGLFD

image



Background

Alternative Names	60S acidic ribosomal protein P2 (Large ribosomal subunit protein P2) (Renal carcinoma antigen NY-REN-44)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use