


Product Name: Recombinant Human TPM1 Protein
Catalog #: PEH6913



Summary

Name	TPM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09493
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDAIKKKMQM LKLDKENALD RAEQAEADKK AAEDRSKQLE DELVSLQKKL KGTEDELDKY SEALKDAQEK LELAEEKATD AEADVASLNR RIQLVEEELD RAQERLATAL QKLEEAEEKAA DESERGMKVI ESRAQKDEEK MEIQEIQLKE AKHIAEDADR KYEEVARKLV IIESDLERAE ERAELSEGKC AELEELKTV TNNLKSLEAQ AEKYSQKEDR YEEEIKVLSD KLKEAETRAE FAERSVTKLE KSIDDLEDEL YAQKLKYKAI SEELDHALND MTSI

image



Background

Alternative Names	Tropomyosin alpha-1 chain (Alpha-tropomyosin) (Tropomyosin-1)
--------------------------	---

Product Name: Recombinant Human TPM1 Protein
Catalog #: PEH6913



Note

For Research Use Only , Not for Diagnostic Use