


Product Name: Recombinant Human HNRNPA1 Protein
Catalog #: PEH6914



Summary

Name	ROA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09651
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSKSESPKEP EQLRKFIGG LSFETTESL RSHFEQWGTL TDCVVMRDPN TKRSRGFGFV TYATVEEVDA AMNARPHKVD GRVVEPKRAV SREDSQRPGA HLTVKKIFVG GIKEDTEHH LRDYFEQYGK IEVIEIMTDR GSGKKRGFAF VTFDDHDSVD KIVIQKYHTV NGHNCVRKA LSKQEMASAS SSQRGRSGSG NFGGGRGGG GGNDFNFRGG NFSGRGGFGG SRGGGGYGGG GDGYNFGND GGYGGGGPGY SGGSRGYGSG GQGYGNQSG YGGSGSYDSY NNGGGGGFGG GSGSNFGGGG SYNDFGNYN QSSNFGPMKG GNFGGRSSGP YGGGGQYFAK PRNQGGYGGG SSSSYGSGR RF

image



Background

Alternative Names	Heterogeneous nuclear ribonucleoprotein A1 (hnRNP A1) (Helix-destabilizing
--------------------------	--

Product Name: Recombinant Human HNRNPA1 Protein
Catalog #: PEH6914



protein) (Single-strand RNA-binding protein) (hnRNP core protein A1)
[Cleaved into: Heterogeneous nuclear ribonucleoprotein A1, N-terminally processed]

Note

For Research Use Only , Not for Diagnostic Use