


Product Name: Recombinant Human PSMB1 Protein
Catalog #: PEH6916



Summary

Name	PSB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20618
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RFSPYVFNGG TILAIAGEDF AIVASDTRLS EGFSIHDRDS PKCYKLTDKT VIGCSGFHGD CLTLTKIEA RLKMYKHSNN KAMTTGAIAA MLSTILYSRR FFPYVYVNI GGLDEEGKGA VYSFDPVGSY QRDSFKAGGS ASAMLQPLLD NQVGFKNMQN VEHVPLSLDR AMRLVKDVFI SAAERDVYTG DALRICIVTK EGIREETVSL RKD

image



Background

Alternative Names	Proteasome subunit beta type-1 (EC 3.4.25.1) (Macropain subunit C5) (Multicatalytic endopeptidase complex subunit C5) (Proteasome component C5) (Proteasome gamma chain)
--------------------------	--

Product Name: Recombinant Human PSMB1 Protein
Catalog #: PEH6916



Note

For Research Use Only , Not for Diagnostic Use