


Product Name: Recombinant Human MPST Protein
Catalog #: PEH6918



Summary

Name	THTM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P25325
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASPQLCRAL VSAQWVAEAL RAPRAGQPLQ LLDASWYLPK LGRDARREFE ERHIPGAFF DIDQCSDRTS PYDHMLPGAE HFAEYAGRLG VGAATHVVIY DASDQGLYSA PRVWWMFRAF GHHAVSLLDG GLRHWLRQNL PLSSGKSQPA PAEFRAQLDP AFIKTYEDIK ENLESRRFQV VDSRATGRFR GTEPEPRDGI EPGHIPGTVN IPFTDFLSQE GLEKSPEEIR HLFQEKKVDL SKPLVATCGS GVTACHVALG AYLCGKPDVP IYDGSWVEWY MRARPEDVIS EGRGKTH

image



Background

Alternative Names 3-mercaptopyruvate sulfurtransferase (MST) (EC 2.8.1.2)

Product Name: Recombinant Human MPST Protein
Catalog #: PEH6918



Note

For Research Use Only , Not for Diagnostic Use