


Product Name: Recombinant Human TMOD1 Protein
Catalog #: PEH6921



Summary

Name	TMOD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P28289
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RDLDEDEILG ALTEELRTL ENELDELPD NALLPAGLRQ KDQTTKAPTG PFKREELLDH LEQAKEFKD REDLVPYTGE KRGKVWVPKQ KPLDPVLESV TLEPELEEAL ANASDAELCD IAAILGMHTL MSNQQYYQAL SSSSIMNKEG LNSVIKPTQY KPVPDEEPNS TDVEETLERI KNNDPKLEEV NLNNIRNIPI PTLKAYAEAL KENSIVKKFS IVGTRSNDPV AYALAEMLKE NKVLKTLNVE SNFISGAGIL RLVEALPYNT SLVEMKIDNQ SQPLGNKVEM EIVSMLEKNA TLLKFGYHFT QQGPRLRASN

image



Background

Alternative Names	Tropomodulin-1 (Erythrocyte tropomodulin) (E-Tmod)
--------------------------	--

Product Name: Recombinant Human TMOD1 Protein
Catalog #: PEH6921



Note

For Research Use Only , Not for Diagnostic Use