

Product Name: Recombinant Human HIBADH Protein
Catalog #: PEH6922



Summary

Name	3HIDH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P31937
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASKTPVGFIG LGNMGNPMAK NLMKHGYPLI IYDVFPDACK EFQDAGEQVV SSPADVAEKA DRIITMLPTS INAIEAYSGA NGILKKVKKG SLLIDSSTID PAVSKELAKE VEKMGAVFMD APVSGGVGAA RSGNLTFMVG GVEDEF AAAQ ELGCMGSNV VYCGAVGTGQ AAKICNNMLL AISMIGTAEA MNLGIRLGLD PKLLAKILNM SSGRCWSSDT YNPVPGVMDG VPSANNYQGG FGTTLMAKDL GLAQDSATST KSPILLGSLA HQIYRMMCAK GYSKKDFSSV FQFLREEETF

image



Background

Alternative Names 3-hydroxyisobutyrate dehydrogenase, mitochondrial (HIBADH) (EC 1.1.1.31)

Product Name: Recombinant Human HIBADH Protein
Catalog #: PEH6922



Note

For Research Use Only , Not for Diagnostic Use