


Product Name: Recombinant Human NNMT Protein
Catalog #: PEH6927



Summary

Name	NNMT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40261
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MESGFTSKDT YLSHFNPRDY LEKYYKFGSR HSAESQILKH LLKNLFKIFC LDGVKGDLLI DIGSGPTIYQ LLSACESFKE IVVTDYSDQN LQELEKWLKK EPEAFDWSPV VTYVCDLEGN RVKGPEKEEK LRQAVKQVLK CDVTQSQPLG AVPLPPADCV LSTLCLDAAC PDLPTYCRAL RNLGSLLKPG GFLVIMDALK SSYYMIGEQQ FSSLPLGRE AVEAAVKEAGY TIEWFEVISQ SYSSTMANNE GLFSLVARKL SRPL

image



Background

Alternative Names	Nicotinamide N-methyltransferase (EC 2.1.1.1)
--------------------------	---

Product Name: Recombinant Human NNMT Protein
Catalog #: PEH6927



Note

For Research Use Only , Not for Diagnostic Use