


Product Name: Recombinant Human TPI1 Protein
Catalog #: PEH6936



Summary

Name	TPIS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P60174-3
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDGEAEFHF AALYISGQWP RLRADTDLQR LGSSAMAPSR KFFVGGNWKM NGRKQSLGEL IGTLNAAKVP ADTEVVCAPP TAYIDFARQK LDPKIAVAAQ NCYKVTNGAF TGEISPGMIK DCGATWVVLG HSERRHVFGE SDELIGQKVA HALAEGLGVI ACIGEKLDER EAGITEKVVF EQTKVIADNV KDWSKVVLAY EPVWAIGTGK TATPQQAQEV HEKLRGWLKS NVSDAVAQST RIIYGGSVTG ATCKELASQP DVD

image



Background

Alternative Names	Triosephosphate isomerase (TIM) (EC 5.3.1.1) (Methylglyoxal synthase) (EC 4.2.3.3) (Triose-phosphate isomerase)
--------------------------	---

Product Name: Recombinant Human TPI1 Protein
Catalog #: PEH6936



Note

For Research Use Only , Not for Diagnostic Use