


Summary

Name	ELAV1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15717
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSNGYEDHMA EDCRGDIGRT NLIVNYLPQN MTQDELRS LF SSIGEVESAK LIRDKVAGHS LYG FVNYVT AKDAERAI NT LNGLRLQSKT IKVSYARPSS EVIKDANLYI SGLPRTMTQK DVEDMFSRFG RIINSRVLVD QTTGLSRGVA FIRFDKRSEA EEAITSFN GH KPPGSSEPIT VKFAANPNQN KNVALLS QLY HSPARRFGGP VHHQAQRFRF SPMGVDHMSG LSGVNVPGNA SSGWCIFIYN LGQDADEGIL WQMFGPFGAV TNVKVIRDFN TNKCKGFGFV TMTNYEEAAM AIASLNGYRL GDKILQVSFK TNKSHK

image



Background

Alternative Names ELAV-like protein 1 (Hu-antigen R) (HuR)

Product Name: Recombinant Human ELAVL1 Protein
Catalog #: PEH6947



Note

For Research Use Only , Not for Diagnostic Use