

Product Name: Recombinant Human HOGA1 Protein
Catalog #: PEH6948



Summary

Name	HOGA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XE5
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASGEGKKVDI AGIYPPVTTT FTATAEVDYG KLEENLHKLK TFPFRGFVVQ GSDGEPFLT SSERLEVSR VRQAMPKNRL LLAGSGCEST QATVENTVSM AQVGADAAMV VTPCYRGRM SSAALIIHHT KVADLSPIPV VLYSVPANTG LDLPVDAVVT LSQHPNIVGM KDSGGDVTRI GLIVHKTRKQ DFQVLAGSAG FLMASALGA VGGVCALANV LGAQVCQLER LCCTGQWEDA QKLQHRLIEP NAAVTRRFGI PGLKKIMDWF GYYGGPCRAP LQELSPAEEE ALRMDFTSNG WL

image



Background

Alternative Names 4-hydroxy-2-oxoglutarate aldolase, mitochondrial (EC 4.1.3.16) (Dihydrodipicolinate synthase-like) (DHDPS-like protein) (Probable 2-keto-4-hydroxyglutarate aldolase) (Probable KHG-aldolase) (Protein 569272)

Product Name: Recombinant Human HOGA1 Protein
Catalog #: PEH6948



Note

For Research Use Only , Not for Diagnostic Use