


Summary

Name	XXLT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NBI6
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATKRLKEARA GAPAAPSPPA LELARGSVAP APGAKAKSLE GGGAGPVDYH LLMMFTKAEH NAALQAKARV ALRSLRLAK FEAEVNLNH FVSEEASREV AKGLLRELLP PAAGFKCKVI FHDVAVLTDK LFPIVEAMQK HFSAGLGTY SDSIFFLSVA MHQIMPKEIL QIIQLDLK FKTNIRELFE EFDSFLPGAI IGIAREMQPV YRHTFWQFRH ENPQTRVGGP PPEGLPGFNS GVMLLNLEAM RQSPLYRLL EPAQVQQLAD KYHFRGHLGD QDFFTMIGME HPKLFHVLDC TWNRLCTWW

image



Background

Alternative Names	Xyloside xylosyltransferase 1 (EC 2.4.2.62) (UDP-xylose:alpha-xyloside alpha-1,3-xylosyltransferase)
--------------------------	--

Product Name: Recombinant Human XXYLT1 Protein
Catalog #: PEH6949



Note

For Research Use Only , Not for Diagnostic Use