


Product Name: Recombinant Human NEK7 Protein
Catalog #: PEH6950



Summary

Name	NEK7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDX7
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDEQSQGMQG PPVPQFQPQK ALRPDMGYNT LANFRIEKKI GRGQFSEVYR AACLLDGVPV ALKKVQIFDL MDAKARADCI KEIDLLKQLN HPNVIKYYAS FIEDNELNIV LELADAGDLS RMIKHFKKQK RLIPERTVWK YFVQLCSALE HMHSRRVMHR DIKPANVFIT ATGVVKLGDL GLGRFFSSKT TAAHSLVGTP YYMSPERIHE NGYNFKSDIW SLGCLLYEMA ALQSPFYGDK MNLYSLCKKI EQCDYPPLPS DHYSEELRQL VNMCPNPDP KRPDVTYVYD VAKRMHACTA SS

image



Background

Alternative Names	Serine/threonine-protein kinase Nek7 (EC 2.7.11.1) (Never in mitosis A-related kinase 7) (NimA-related protein kinase 7)
--------------------------	--

Product Name: Recombinant Human NEK7 Protein
Catalog #: PEH6950



Note

For Research Use Only , Not for Diagnostic Use